

# **Northwest**

**Life Science Specialties, LLC**

**Premier Products for Superior Life Science Research**

***NWLSSTM***

***Nitric Oxide (Nitrate/Nitrite)***

***Enzymatic Assay***

---

**Product NWK-EN001**

***For Research Use Only***

Enzymatic assay system for measurement of nitric oxide in biological samples. This method employs a nitrate reductase to reduce nitrates to nitrate before measuring using Griess reagent.

## Table of Contents

---

<b>Section</b>	<b>Page</b>
Introduction	3
Intended Use	3
Test Principle	3
General Specifications	3
Kit Contents	4
Required Materials Not Provided	4
Required Instrumentation	4
Warnings, Limitations, Precautions	4
Storage Instructions	4
Assay Preparation	4
Reagent Preparation	5
Standard Preparation	5
Sample Handling/Preparation	6
Assay Protocol	6
Data Analysis	7
Performance Details	8
References	9
Statement of Limited Warranty	9
User Notes	10

**Introduction:**

Nitric Oxide (NO) is a biologically relevant free radical and cell signaling molecule. In normal cellular physiology it was originally characterized as an endothelial relaxation factor (EDRF) and is most commonly associated with its role in regulation of vasodilation. From a pathogenic standpoint the free radical nitric oxide (NO•) is known to combine with superoxide (O<sub>2</sub>•) to form peroxynitrite (ONOO<sup>-</sup>) which is a potent oxidant. Since ONOO<sup>-</sup> (NO<sub>3</sub><sup>-</sup>) is essentially an unstable isomer of nitrate, it can react with tyrosine residues in proteins to create nitrotyrosine, a biomarker of oxidative mediated modification of proteins.



Nitric oxide degrades rapidly to nitrate (NO<sub>3</sub>) and nitrite (NO<sub>2</sub>) in aqueous biological systems. When NO<sub>3</sub> is reduced to NO<sub>2</sub> measurement of nitrite provides a simple, indirect means to quantify nitric oxide produced in a wide array of experimental model systems.

**Intended Use:**


This NWLSS™ Nitric Oxide (Nitrate/Nitrite) Enzymatic Assay kit is designed for use in quantifying nitric oxide (NO) in biological samples.

**Test Principle:**

The test method is based on measurement of total sample nitrites using Griess Reagent. Since Griess reagent does not detect nitrate NO<sub>3</sub>, This kit employs nitrate reductase for reduction of nitrate to nitrite (Fig. 1).

**Figure 1: Nitrate Reduction****Figure 2: Sulfanilamide reaction with HNO<sub>2</sub>**

After the appropriate sample reduction step nitrites are converted to nitrous acid (HNO<sub>2</sub>) in Sulfanilamide, HCl solution. HNO<sub>2</sub> subsequently reacts with sulfanilamide to form sulfanilamide diazonium salt (Fig 2).

**Figure 3: Diazonium Salt Reaction with NED**

Diazonium salt is then reacted with N-(1-Naphthyl)-ethylenediamine (NED) to produce a chromophore directly measurable at 540 nm (Fig. 3).

**General Specifications:**

Format: 96 wells

Number of tests: Triplicate = 24  
Duplicate = 40

Specificity: Nitric Oxide as total Nitrite

Sensitivity: 1 pmol/mL or 1 $\mu$ M in the assay

Sample Volume: 10 - 50  $\mu$ M

Effective Range: 1  $\mu$ M - 100  $\mu$ M

**Kit Contents**

Microplate (96 well clear, low binding, flatbottom)	1 X 96 wells
MOPS Buffer (Sample Dilution Buffer)	1 X 25 mL
Nitrate Standard (500 $\mu$ M KNO <sub>3</sub> )	1 X 1.5 mL
Nitrate Reductase (Lyophilized)	1 Unit
Nitrate Reductase Buffer	1.5 mL
NADH (Lyophilized)	2 mg
Reagent A:	1 X 7 mL
(Sulfanilamide (p-Aminobenzenesulfonamide) in 3N HCl)	
Reagent B:	1 X 7 mL
(N-(1-Naphthyl) ethylenediamine dihydrochloride in deionized H <sub>2</sub> O)	

**Required Materials Not Provided:**

Adjustable pipettes with disposable tips (range of 10  $\mu$ L to 1,000  $\mu$ L)  
Deionized water.

Microcentrifuge tubes.

Vortex Mixer.

**Required Instrumentation:**

Microtiter plate reader with 540 nm capability.

**Warnings, Limitations, Precautions:**

Individual components may be harmful if swallowed, inhaled or absorbed through the skin. Contact should be minimized through the use of gloves and standard good laboratory practices. If contact with skin or eyes occurs, rinse the site immediately with water and consult a physician.

**Storage Instructions:**

Upon receipt, store Nitrate Reductase enzyme at -20 °C. Store NADH in the dark at room temperature. Store all other components at 4 °C until immediately before use.

**Assay Preparation**

1. Determine the number of wells required to assay standards, samples and controls for the appropriate replicates.
2. Create an assay template showing positioning of standards, controls and samples. Include blank wells also.

**Reagent Preparation:*****Nitrate Reductase***

Reconstitute with 1.0 mL *Nitrate Reductase Buffer* and incubate at room temperature for 20 minutes. Vortex gently at 0, 10, and 20 minutes.

***NADH***

Add 1.28 mL deionized water to the vial to create 2 mM NADH Reagent. Label as *Working NADH Reagent*.

All other reagents are supplied ready to use.

**Standard Preparation:**

1. Standard Supplied: 1.5 mL of 500  $\mu\text{M}$   $\text{KNO}_3$
2. Label 8 tubes ( $S_7$  to  $S_0$ ) as 50, 25, 12.5, 6.25, 3.13, 1.56, 0.78 and 0  $\mu\text{M}$ .
3. In tube 8 ( $S_7$ ): Add 100  $\mu\text{L}$  of 500  $\mu\text{M}$   $\text{NaNO}_2$  Standard as supplied to 900  $\mu\text{L}$  deionized  $\text{H}_2\text{O}$ . Mix well.  $S_7$  is now 50  $\mu\text{M}$ .
4. Perform a serial dilution by transferring 400  $\mu\text{L}$  of 50  $\mu\text{M}$   $S_7$  to  $S_6$  then 400  $\mu\text{L}$   $S_6$  to  $S_5$  and so on through to  $S_1$  to create all standards down to 0.78  $\mu\text{M}$ . Leave  $S_0$  as deionized  $\text{H}_2\text{O}$  only zero or blank.

**Sample Handling/Preparation:***Samples with low levels of nitrite:*

If the nitrite in a sample is too low, sample volume may be increased while decreasing the volume of buffer added.

*Samples with high levels of protein:*

Excess protein in samples can produce a precipitate that may interfere with the measurement of NO. If there is visible precipitate after addition of Reagent A, if possible, samples can be diluted to reduce the amount of protein present and retested.


**Note:** If sample proteins cannot be diluted enough to allow for accurate testing using this assay, we recommend using our Non-Enzymatic Nitric Oxide Assay, product number NWK-NN001.

**Assay Procedure:**

1. Add 85  $\mu$ L of each **standard** and **sample** to replicate wells according to assay template.
2. Add 10  $\mu$ L reconstituted *Nitrate Reductase* to each well.
3. Add 10  $\mu$ L of *Working NADH Reagent* to each well and shake the plate for 20 minutes at room temperature.
4. Add 50  $\mu$ L of **Reagent A** to each well and shake briefly.
3. Add 50  $\mu$ L of **Reagent B** to each well and shake for 5 minutes at room temperature.
5. Read and record absorbance at 540 nm.

**Data Analysis**

1. Average the  $A_{540}$  values for each replicate of sample, standard and blank.
2. Subtract the average  $A_{540}$  value of the blank wells from the average  $A_{540}$  for each standard and sample replicate.
3. Plot a standard curve as concentration  $\text{NO}_2$  ( $\mu\text{M}$  Nitric Oxide Equivalents) vs. Absorbance ( $A_{540}$ ).



4. Calculate the concentration of nitrite in each sample using the equation as derived from the standard curve generated.

Notes: Each plate tested must have its own standard curve. The Example standard curve above is for illustration purposes only.

Values obtained must be multiplied by dilution factors incurred during sample precipitation and reduction steps and any other dilutions that may have been necessary.

Since  $1 \mu\text{M}$  Nitric Oxide Equivalents =  $1 \text{ pmol/mL}$   $\text{NO}_3$ , the standard curve can also be plotted as  $\text{NO}_3$  concentration ( $\text{pmol/mL}$ ) vs. Absorbance ( $A_{540}$ ).

**REFERENCES**

1. Schmidt, H.H., et. al., (1995) *Biochemica* 2:22-23
2. Campbell, E. R., et. Al., (2000) *American Laboratory* February, 90-92.

**Statement of Limited Warranty:**


Northwest Life Science Specialties, LLC (NWLSS) makes no guarantee of any kind, expressed or implied, that extends beyond the description of the material in this kit, except that they will meet our specifications at the time of delivery. Customer's remedy and NWLSS' sole liability is limited to, at NWLSS' option, refund of the purchase price, or the replacement of material not meeting our specification. By acceptance of our product, customer assumes all liability and will indemnify and hold NWLSS harmless for the consequence of this product's use or misuse by the customer, its employees, or others. Refund or replacement is conditioned of customer notifying NWLSS within twenty-one (21) days of the receipt of product. Failure to give notice within 21 days shall constitute a waiver by the customer of all claims hereunder with respect to said product.





User Notes:

User Notes:



# **Northwest**

**Life Science Specialties, LLC**

**5131 NE 94th Avenue, Suite 101**

**Vancouver, WA 98662**

**Phone 360-449-3091 or Toll Free: 888-449-3091**

**Fax 360-449-3092**

**E-mail: [sales@nwlifescience.com](mailto:sales@nwlifescience.com)**