


Northwest

Life Science Specialties, LLC

Anti Myeloperoxidase monoclonal Clone: 1A1

Catalog Number:	NWA-MPOM01A1
Description:	Murine IgG _{1κ} anti Hu-Myeloperoxidase 100 µg (100 µL at 1mg/mL)
Specificity:	Specific for human myeloperoxidase
Format:	100 uL IgG in phosphate buffered saline (PBS) and 50% glycerol
Antigen:	Human Myeloperoxidase purified from sputum.
Suggested Use:	ELISA Immunoprecipitation (2 µL/400 µL Lysate) <i>User should determine optimum titer for each application.</i>
Cross Reactivity:	Human
Storage:	Store at -20 °C upon arrival.
Stability:	Store for 1 year at -20 degree from date of shipment



Immunoprecipitation Analysis of HL-60 Cell Lysates. Lane 1: Input, Lane 2: Precipitates

- References:**
- 1) Winterbourn CC, et al, Curr Opin Hematol. 200; vol.7(1): pp.53-8.
 - 2) Hansson M, et al, Arch Biochem Biophys. 2006; vol.445(2): pp.214-24.
 - 3) Nicholls SJ, Hazen SL. Arterioscler Thromb Vasc Biol. 2005; vol.25(6): pp.1102-11.
 - 4) Klebanoff SJ. J Leukoc Biol. 2005; vol.77(5): pp.598-625.

For Research Use Only

Northwest Life Science Specialties, LLC
16420 S McGillivray, Suite 103, PMB 106, Vancouver, WA 98683
Phone: (360) 449-3091, (888) 449-3091, Fax: (360) 449-3092, E-mail: info@nwlifescience.com