


Northwest

Life Science Specialties, LLC

Anti Peroxiredoxin I Monoclonal Clone 2A4

Catalog Number:	NWA-PRX102A4
Description:	Murine IgG _{1κ} anti Hu-Peroxiredoxin I (100 µL per Vial)
Specificity:	Specific for human Peroxiredoxin I
Format:	100 uL IgG in HEPES with 0.15M NaCl, 0.01% BSA, 0.03% sodium azide, and 50% glycerol
Antigen:	Recombinant Human Peroxiredoxin I from E. coli.
Suggested Use:	IP (1-2 µL/400uL Lysate) <i>User should determine optimum titer for each application.</i>
Cross Reactivity:	Human & Rat
Storage:	Store at -20 °C upon arrival.
Stability:	Store for 1 year at -20 degree from date of shipment



IMMUNOPRECIPITATION ANALYSIS of HeLa cell lysates:
Lane 1: Input Lane 2: Precipitates

- References:**
- (1) Wood, Z. A. et al. (2003) Trends Biochem Sci. 28(1):32-40.
 - (2) Rhee Sue Goo, et al (2001) IUBMB life 52:35-41
 - (3) Min Hee Choi, et al (2005) Nature letters 435(19) :347-353

For Research Use Only

Northwest Life Science Specialties, LLC
16420 S McGillivray, Suite 103, PMB 106, Vancouver, WA 98683
Phone: (360) 449-3091, (888) 449-3091, Fax: (360) 449-3092, E-mail: info@nwlifescience.com