


Northwest

Life Science Specialties, LLC

Anti Peroxiredoxin III Monoclonal Clone 2B11

Catalog Number:	NWA-PRX32B11
Description:	Murine IgG _{2ak} anti Hu-Peroxiredoxin III (100 µL per Vial)
Specificity:	Specific for human Peroxiredoxin III
Format:	100 uL IgG in HEPES with 0.15M NaCl, 0.01% BSA, 0.03% sodium azide, and 50% glycerol
Antigen:	Recombinant Human Peroxiredoxin III from E. coli.
Suggested Use:	ELISA IP (1-2 µL/400uL Lysate) <i>User should determine optimum titer for each application.</i>
Cross Reactivity:	Human & Rat
Storage:	Store at -20 °C upon arrival.
Stability:	Store for 1 year at -20 degree from date of shipment



IMMUNOPRECIPITATION ANALYSIS of 293T cell lysates:
Lane 1: Input Lane 2: Precipitates

References: 1) Wood, Z. A. et al. (2003) Trends Biochem Sci. 28(1):32-40.

For Research Use Only

Northwest Life Science Specialties, LLC
16420 S McGillivray, Suite 103, PMB 106, Vancouver, WA 98683
Phone: (360) 449-3091, (888) 449-3091, Fax: (360) 449-3092, E-mail: info@nwlifescience.com